

Pen Needles & Insulin Syringes

PREMIUM QUALITY NEEDLES FOR DIABETES CARE



PEN NEEDLE
UltiGuard® SafePack

UltiCare

You play a critical role in your community's health.

Pharmacists like you play critical roles in the care continuum. After all, many of your patients see you more often than they see their doctors.

When you select UltiGuard Safe Pack, you give your patients more for the same copay and bring in higher revenue and margins for your business. Further, as most of your patients don't buy sharps containers, dispensing UltiGuard Safe Pack means you do your part for sharps stewardship in your community.

With UltiGuard Safe Pack, you can do the right thing for your patients, your community, and your pharmacy.

**It's Your Choice
Choose UltiCare**



Why UltiGuard Safe Pack

When you choose the all-in-one UltiGuard Safe Pack, you give your patients premium, 100-count pen needles and a sharps container, making sharps disposal convenient & affordable.

With the FDA-cleared UltiGuard Safe Pack system, your patients don't need to buy a separate sharps container—and the UltiGuard Safe Pack has the same copay as other options. With UltiGuard Safe Pack, you give your patients more for no additional cost.



Not shown to actual size. 4.5" x 6"

How UltiGuard Safe Pack Works

UltiGuard Safe Pack has an upper and a lower chamber. Patients get new pen needles from the lower chamber and safely discard used pen needles in the upper chamber. Pen Needle UltiGuard Safe Pack safeguards used needles once patients turn the wheel in the container's side.

When patients finish their needle supplies, they can dispose the UltiGuard Safe Pack according to their state or local guidelines.

UltiCare® Pen Needles

Available in the UltiGuard Safe Pack — the only dual-purpose system that allows you to dispense a new needle and safely dispose a used needle. Our premium quality needles offer all the comfort features of the leading brand with our patented all-in-one solution specially designed for in-home use.

The Safe Pack contains 100 needles and is available in Micro (4mm), Mini (5mm), Mini (6mm), Short (8mm), and Original (12.7mm) sizes. Also available in 30, 50, and 100 count boxes. UltiCare pen needles fit all pen injectors sold in the U.S.



ULTIGUARD 100-COUNT	BOX 30-COUNT	BOX 50-COUNT	BOX 100-COUNT
Micro 4mm (5/32") x 32G			
NDC#08222.0954.35 ABC# 10103231	NDC#08222.7100.79 ABC# 10162944	NDC#08222.0954.59 ABC# 10103272	NDC#08222.7100.31 ABC# 10210612
Mini 5mm (3/16") x 31G			
NDC#08222.0955.34 ABC# 10213659			NDC#08222.7101.30 ABC# 10213366
Mini 6mm (1/4") x 32G			
NDC#08222.0956.26 ABC# 10213658			NDC#08222.7101.47 ABC# 10213657
Mini 6mm (1/4") x 31G			
NDC#08222.0956.33 ABC# 10051644	NDC#08222.7100.86 ABC# 10162943	NDC#08222.0956.57 ABC# 10095525	NDC#08222.7100.48 ABC# 10209198
Short 8mm (5/16") x 31G			
NDC#08222.0958.31 ABC# 10052938	NDC#08222.7100.93 ABC# 10162942	NDC#08222.0958.55 ABC# 10095526	NDC#08222.7100.55 ABC# 10210615
Original 12.7mm (1/2") x 29G			
NDC#08222.0951.21 ABC# 10024752			NDC#08222.7100.62

UltiCare Tuberculin/Flu Safety Syringes

The UltiCare Tuberculin/Flu Safety Syringe, with a permanently attached needle, offers an audible confirmation of the lock position.

1mL	25.4mm (1") x 25G	1mL	15.9mm (5/8") x 25G
NDC#08222.2511.07 ABC# 10098795		NDC#08222.2515.89	

UltiCare Sterile Alcohol Swabs

The ultra-soft and thick UltiCare Sterile Alcohol Swabs optimally clean the skin prior to injection to prevent infection.

BOX 100-COUNT	BOX 200-COUNT
NDC#08222.4400.06 ABC# 10152307	NDC#08222.4400.13 ABC# 10152305

UltiCare® Insulin Syringes

UltiCare insulin syringes, similar to our pen needles, are as comfortable as the leading brand and offer the same UltiGuard Safe Pack dual-purpose system. The Safe Pack contains 100 syringes and is available in 3/10mL, 1/2mL, and 1mL sizes.

Also available in a 100-count box.



3/10mL (30-UNIT)

ULTIGUARD 100-COUNT	BOX 100-COUNT
6mm (1/4") x 31g	
NDC#08222.9101.03	NDC#08222.9100.66 ABC# 10164469
6mm (1/4") x 31g with half-unit scale	
NDC#08222.9101.10	NDC#08222.9100.73 ABC# 10164468
8mm (5/16") x 31g	
NDC#08222.0743.93 ABC# 10018365	NDC#08222.0943.91 ABC# 10030524
8mm (5/16") x 30g	
	NDC#08222.0933.94 ABC# 10030561
12.7mm (1/2") x 30g	
NDC#08222.0733.58 ABC# 10018359	NDC#08222.0933.56 ABC# 10004811
12.7mm (1/2") x 29g	
	NDC#08222.0923.97 ABC# 10030863

1/2mL (50-UNIT)

ULTIGUARD 100-COUNT	BOX 100-COUNT
6mm (1/4") x 31g	
NDC#08222.9101.27	NDC#08222.9100.80 ABC# 10164464
8mm (5/16") x 31g	
NDC#08222.0745.91 ABC# 10030851	NDC#08222.0945.99 ABC# 10034670
8mm (5/16") x 30g	
	NDC#08222.0935.92 ABC# 10030562
12.7mm (1/2") x 30g	
NDC#08222.0735.56 ABC# 10018354	NDC#08222.0935.54 ABC# 10007636
12.7mm (1/2") x 29g	
	NDC#08222.0925.95 ABC# 10030860
12.7mm (1/2") x 28g	
	NDC#08222.0825.89

1mL (100-UNIT)

ULTIGUARD 100-COUNT	BOX 100-COUNT
6mm (1/4") x 31g	
NDC#08222.9101.34	NDC#08222.9100.97 ABC# 10164480
8mm (5/16") x 31g	
NDC#08222.0741.95 ABC# 10030852	NDC#08222.0941.93 ABC# 10030554
8mm (5/16") x 30g	
	NDC#08222.0931.96 ABC# 10030568
12.7mm (1/2") x 30g	
NDC#08222.0731.50 ABC# 10018364	NDC#08222.0931.58 ABC# 10025222
12.7mm (1/2") x 29g	
	NDC#08222.0921.99 ABC# 10034664
12.7mm (1/2") x 28g	
	NDC#08222.0821.83 ABC# 10003989

UltiMed Incorporated

is a privately held company that manufactures premium quality insulin syringes, pen needles, and safety syringes for the diabetes and animal health markets. The company was founded in 1988 to create high quality products at a better price.

The UltiGuard Safe Pack by UltiMed, is the only FDA cleared dual-purpose system that combines premium quality needles and a sharps container. UltiMed's exclusive all-in-one solution makes the company the world leader in safe needle disposal.

UltiMed is based in Minneapolis, Minnesota with manufacturing facilities in the United States and abroad.



ULTIMED
INCORPORATED

ulticare.com

1.877.854.3434