

DID YOU KNOW?

Home-generated sharps disposal regulations have changed in California effective January 7th, 2022.

Under California Senate Bill 212, pharmacies must provide or initiate distribution of a sharps waste container and mail-back materials.

(SB 212, Jackson. Solid waste: pharmaceutical and sharps waste stewardship)

The UltiGuard Safe Pack Sharps Container & Mail-Back Disposal Kit is the only fully compliant all-in-one solution.



**COMPLETE
3-IN-1
SOLUTION**



1

100 pen needles
or syringes



2

FDA-cleared
sharps container



3

Prepaid return
shipping box and
mail-back materials

GIVE YOUR PATIENTS MORE



FOR THE SAME COPAY AS PEN NEEDLES
OR SYRINGES ALONE



WITH NO CHANGE IN PHARMACY WORKFLOW



IN ONE COMPLETE, PREPAID MAIL-BACK KIT

UltiGuard® SafePack
Sharps Container & Mail-Back Disposal Kit

See how we make it easy to stay compliant:
www.disposemyneedles.com

How It Works

1. Box

Place the Safe Pack into the reclosable plastic bag, seal it, then place the sealed bag into the prepaid return shipping box.

3. Sign

Remove the tracking manifest from the plastic bag on the side of the shipping box and fill it in completely. Retain the last page for your records, then place the completed manifest back into the bag.

2. Seal

Close the brown inner flap. Then close the outer white shipping box following the numbers. Seal the closure tab with the provided tape strip.

4. Send

Take to the nearest U.S. Postal Service Office or place in your mailbox for pickup.



Pen Needle & Syringe Sharps Container & Mail-Back Disposal Kit Products

Available in all standard pen needle and insulin syringe sizes. Each kit includes 100 pen needles or syringes, an FDA-cleared sharps container and a prepaid return shipping box. The kit can be ordered from all national wholesalers including AmerisourceBergen, Cardinal Health, and McKesson.



4mm x 32G Micro
NDC 08222-1254-39



5mm x 31G Mini
NDC 08222-1255-38



6mm x 32G Mini
NDC 08222-1256-20



8mm x 31G Short
NDC 08222-1258-35



12.7mm x 29G Original
NDC 08222-1251-25



**3/10 mL/cc 8mm
(5/16'') x 31G**
NDC 08222-1243-95



**1/2 mL/cc 8mm
(5/16'') x 31G**
NDC 08222-1245-93



**1 mL/cc 8mm
(5/16'') x 31G**
NDC 08222-1241-97



**3/10 mL/cc 12.7mm
(1/2'') x 30G**
NDC 08222-1233-50



**1/2 mL/cc 12.7mm
(1/2'') x 30G**
NDC 08222-1235-58



**1 mL/cc 12.7mm
(1/2'') x 30G**
NDC 08222-1231-52