

Upgrade Your Syringe.

UltiGuard Safe Pack U-100 Insulin Syringes & UltiCare U-100 Insulin Syringes.

UltiMed offers two syringe products — the UltiGuard Safe Pack and UltiCare. Both product lines offer premium syringes, but the UltiGuard Safe Pack also comes with a built-in sharps container.

Get More For the Same Copay.

UltiGuard Safe Pack is the only dual-purpose system that allows you to dispense a new syringe and safely dispose of used syringes from the same container and for the same copay as insulin syringes alone. There's no need to spend more on a separate sharps container. The system includes an FDA-cleared sharps container that meets safe needle disposal guidelines in all 50 states.

Keep Your Home and Community Safe.

Surveys indicate that less than 5% of the more than 7.8 billion sharps devices sold annually in the United States are disposed of in a closed container. Of the remaining 95%, most are deposited unprotected into the household trash, which increases the risk of needle sticks.

With the UltiGuard Safe Pack all-in-one system, you don't have to worry about safe disposal or needle stick injuries. When you choose UltiGuard Safe Pack Insulin Syringes, you choose premium insulin syringes and you choose to keep your home and community safe.

How UltiGuard Safe Pack Works.

UltiGuard Safe Pack has an upper and a lower chamber. New syringes are dispensed from the lower chamber and safely discarded in the upper chamber. UltiGuard Safe Pack secures used syringes once the wheel on the side of the container has been turned.

When you have used all the syringe supply, follow your state or local guidelines for proper disposal.

It's Your Choice.

Your pharmacist chooses the syringe you get with your insulin or other injection medication. If you would like to upgrade your syringes, simply print out page two and three of this document and bring it to your pharmacist.



Product Comparison

PRODUCT	UltiGuard Safe Pack U-100 Insulin Syringes	UltiCare U-100 Insulin Syringes
DESCRIPTION	Premium syringes with a built-in sharps container	Premium syringes
SIZES	Available in all sizes	Available in all sizes
QUANTITY	100	100
BUILT-IN SHARPS CONTAINER	Yes	No
INSURANCE	Covered by most insurance companies	Covered by most insurance companies

It's easy to switch.



UltiMed Insulin Syringe Ordering Information:

UltiMed's two consumer insulin syringe product lines — UltiGuard Safe Pack and UltiCare are available at most pharmacies or can easily be ordered by your pharmacist. Simply print out this page and bring it into your pharmacy to share with your pharmacist.

Name	Needle Length	Gauge	Quantity	NDC	Built-in Sharps Container?
ULTIGUARD SAFE PACK INSULIN SYRINGES					
3/10 mL/cc 6mm x 31G	6mm (1/4")	31G (0.25mm)	100	08222-9101-03	Yes
3/10 mL/cc 6mm x 31G with Half Unit Markings	6mm (1/4")	31G (0.25mm)	100	08222-9101-10	Yes
3/10 mL/cc 12.7mm x 30G	12.7mm (1/2")	30G (0.3mm)	100	08222-0733-58	Yes
3/10 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	08222-0743-93	Yes
1/2 mL/cc 6mm x 31G	6mm (1/4")	31G (0.25mm)	100	08222-9101-27	Yes
1/2 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	08222-0745-91	Yes
1/2 mL/cc 12.7mm x 30G	12.7mm (1/2")	30G (0.3mm)	100	08222-0735-56	Yes
1 mL/cc 6mm x 31G	6mm (1/4")	31G (0.25mm)	100	08222-9101-34	Yes
1 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	08222-0741-95	Yes
1 mL/cc 12.7mm x 30G	12.7mm (1/2")	30G (0.3mm)	100	08222-0731-50	Yes



Alcohol Swabs

The UltiCare two-ply 70% isopropyl alcohol swabs optimally clean the skin prior to injection to prevent infection. Each box contains individually wrapped sterile swabs. Ask your pharmacist today to add this to your order.



UltiMed Insulin Syringe Ordering Information Continued:

UltiMed's two consumer insulin syringe product lines — UltiGuard Safe Pack and UltiCare are available at most pharmacies or can easily be ordered by your pharmacist. Simply print out this page and bring it into your pharmacy to share with your pharmacist.



Name	Needle Length	Gauge	Quantity	NDC	Built-in Sharps Container?
ULTICARE INSULIN SYRINGES					
3/10 mL/cc 6mm x 31G	6mm (1/4")	31G (0.25mm)	100	08222-9100-66	No
3/10 mL/cc 6mm x 31G with Half Unit Marketings	6mm (1/4")	31G (0.25mm)	100	08222-9100-73	No
3/10 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	08222-0943-91	No
3/10 mL/cc 8mm x 30G	8mm (5/16")	30G (0.3mm)	100	08222-0933-94	No
3/10 mL/cc 12.7mm x 30G	12.7mm (1/2")	30G (0.3mm)	100	08222-0933-56	No
3/10 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	08222-0923-97	No
1/2 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	08222-0945-99	No
1/2 mL/cc 6mm x 31G	6mm (1/4")	31G (0.25mm)	100	08222-9100-80	No
1/2 mL/cc 8mm x 30G	8mm (5/16")	30G (0.3mm)	100	08222-0935-92	No
1/2 mL/cc 12.7mm x 30G	12.7mm (1/2")	30G (0.3mm)	100	08222-0935-54	No
1/2 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	08222-0925-95	No
1/2 mL/cc 12.7mm x 28G	12.7mm (1/2")	28G (0.36mm)	100	08222-0825-89	No
1 mL/cc 6mm x 31G	6mm (1/4")	31G (0.25mm)	100	08222-9100-97	No
1 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	08222-0941-93	No
1 mL/cc 8mm x 30G	8mm (5/16")	30G (0.3mm)	100	08222-0931-96	No
1 mL/cc 12.7mm x 30G	12.7mm (1/2")	30G (0.3mm)	100	08222-0931-58	No
1 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	08222-0921-99	No
1 mL/cc 12.7mm x 28G	12.7mm (1/2")	28G (0.36mm)	100	08222-0821-83	No