

Upgrade Your Pen Needle.

UltiGuard Safe Pack Pen Needles & UltiCare Pen Needles.

UltiMed offers two pen needle products — UltiGuard Safe Pack and UltiCare. Both product lines offer premium pen needles, but the UltiGuard Safe Pack also comes with a built-in sharps container.

Get More For the Same Copay.

UltiGuard Safe Pack is the only dual-purpose system that allows you to dispense a new pen needle and safely dispose of used needles from the same container and for the same copay as pen needles alone. There's no need to spend more on a separate sharps container. The system includes an FDA-cleared sharps container that meets safe needle disposal guidelines in all 50 states.

Keep Your Home and Community Safe.

Surveys indicate that less than 5% of the more than 7.8 billion sharps devices sold annually in the United States are disposed of in a closed container. Of the remaining 95%, most are deposited unprotected into the household trash, which increases the risk of needle sticks.

With the UltiGuard Safe Pack all-in-one system, you don't have to worry about safe disposal or needle stick injuries. When you choose UltiGuard Safe Pack Pen Needles, you choose premium pen needles and you choose to keep your home and community safe.

How UltiGuard Safe Pack Works.

UltiGuard Safe Pack has an upper and a lower chamber. New pen needles are dispensed from the lower chamber and safely discarded in the upper chamber. UltiGuard Safe Pack secures used pen needles once the wheel on the side of the container has been turned.

When you have used all the pen needle supply, follow your state or local guidelines for proper disposal.

It's Your Choice.

Your pharmacist chooses the pen needle you get with your insulin or other injection medication. If you would like to upgrade your pen needles, simply print out page two and three of this document and bring it to your pharmacist.





It's easy to switch.

UltiMed Pen Needle Ordering Information:

UltiMed's two consumer pen needle product lines — UltiGuard Safe Pack and UltiCare, are available at most pharmacies and can easily be ordered by your pharmacist. Simply print out this page and bring it into your pharmacy to share with your pharmacist.

Product Comparison

PRODUCT	DESCRIPTION	SIZES	QUANTITY	BUILT-IN SHARPS CONTAINER	INSURANCE	COMPATIBILITY
	UltiGuard Safe Pack Pen Needles	Premium pen needles with a built-in sharps container	Available in all sizes	100	Yes	Covered by most insurance companies. Compatible with most pen injectors devices.
	UltiCare Pen Needles	Premium pen needles	Available in all sizes	30, 50, or 100	No	Covered by most insurance companies. Compatible with most pen injectors devices.

Name	Needle Length	Gauge	Quantity	NDC	Built-in Sharps Container?
ULTIGUARD SAFE PACK PEN NEEDLES					
4mm x 32G Micro	4mm (5/32")	32G (0.23mm)	100	08222-0954-35	Yes
5mm x 31G Mini	5mm (3/16")	31G (0.25mm)	100	08222-0955-34	Yes
6mm x 32G Mini	6mm (1/4")	32G (0.23mm)	100	08222-0956-26	Yes
6mm x 31G Mini	6mm (1/4")	31G (0.25mm)	100	08222-0956-33	Yes
8mm x 31G Short	8mm (5/16")	31G (0.25mm)	100	08222-0958-31	Yes
12.7mm x 29G Original	12.7mm (1/2")	29G (0.33mm)	100	08222-0951-21	Yes



(Continued) UltiCare Pen Needle Ordering Information:

Name	Needle Length	Gauge	Quantity	NDC	Built-in Sharps Container?
ULTICARE PEN NEEDLES					
4mm x 32G Micro					
30 count	4mm (5/32")	32G (0.23mm)	30	08222-7100-79	No
50 count	4mm (5/32")	32G (0.23mm)	50	08222-0954-59	No
100 count	4mm (5/32")	32G (0.23mm)	100	08222-7100-31	No
5mm x 31G Mini					
100 count	5mm (3/16")	31G (0.25mm)	100	08222-7101-30	No
6mm x 32G Mini					
100 count	6mm (1/4")	32G (0.23mm)	100	08222-7101-47	No
6mm x 31G Mini					
30 count	6mm (1/4")	31G (0.25mm)	30	08222-7100-86	No
50 count	6mm (1/4")	31G (0.25mm)	50	08222-0956-57	No
100 count	6mm (1/4")	31G (0.25mm)	100	08222-7100-48	No
8mm x 31G Short					
30 count	8mm (5/16")	31G (0.25mm)	30	08222-7100-93	No
50 count	8mm (5/16")	31G (0.25mm)	50	08222-0958-55	No
100 count	8mm (5/16")	31G (0.25mm)	100	08222-7100-55	No
12.7mm x 29G Original					
100 count	12.7mm (1/2")	29G (0.33mm)	100	08222-7100-62	No



Alcohol Swabs

The UltiCare two-ply 70% isopropyl alcohol swabs optimally clean the skin prior to injection to prevent infection. Each box contains individually wrapped sterile swabs. Ask your pharmacist today to add this to your order.