

Care for your pet. Protect your home.



UltiCare VetRx Syringes & Pen Needles.

UltiMed offers two veterinary insulin syringe product lines — UltiCare VetRx UltiGuard Safe Pack Insulin Syringes and UltiCare VetRx Insulin Syringes. Both product lines offer premium syringes, but the UltiGuard Safe Pack also comes with a built-in sharps container. UltiMed also offers veterinary pen needles: UltiCare VetRx UltiGuard Safe Pack Pen Needles. These pen needles also come with a built-in sharps container.

Upgrade Your Needles & Syringes.

The UltiCare VetRx UltiGuard Safe Pack is the only dual-purpose system that allows you to dispense a new syringe or pen needle and safely dispose of used syringes via the same container. There's no need to spend more on a separate sharps container. The kit is an FDA-cleared sharps container and meets safe needle disposal guidelines in all 50 states.

Keep Your Home and Community Safe.

Surveys indicate that less than 5% of the more than 7.8 billion sharps sold annually in the United States are disposed of in a closed container. Of the remaining 95%, most are deposited unprotected into the household trash, which increases the risk of needle sticks.




With the UltiCare VetRx UltiGuard Safe Pack all-in-one system, you don't have to worry about safe disposal or accidental needle stick injuries. When you choose UltiCare VetRx UltiGuard Safe Pack Insulin Syringes or Pen Needles, you choose premium syringes or pen needles and choose to keep your home and community safe.

How UltiGuard Safe Pack Works.

UltiGuard Safe Pack has an upper and a lower chamber. New syringes are dispensed from the lower chamber and safely discarded in the upper chamber. UltiGuard Safe Pack secures used syringes once the wheel on the side of the container has been turned.

When you have used all the syringe supply, follow your state or local guidelines for proper disposal.

Product Comparison

PRODUCT	DESCRIPTION	SIZES	QUANTITY	BUILT-IN SHARPS CONTAINER	COMPATIBILITY	
	UltiCare VetRx UltiGuard Safe Pack Insulin Syringes	Premium U-40 and U-100 insulin syringes with a built-in sharps container	Available in all sizes	100	Yes	—
	UltiCare VetRx Insulin Syringes	Premium U-40 and U-100 insulin syringes	Available in all sizes	100	No	—
	UltiCare VetRx UltiGuard Safe Pack Pen Needles	Premium pen needles with a built-in sharps container	12.7mm	100	Yes	Compatible with most pen injectors devices.

It's easy to upgrade.

UltiCare VetRx UltiGuard Safe Pack Ordering Information:

UltiMed offers three veterinary UltiGuard Safe Pack product lines — UltiCare VetRx UltiGuard Safe Pack U-100 Insulin Syringes, UltiCare VetRx UltiGuard Safe Pack U-40 Insulin Syringes, and UltiCare VetRx UltiGuard Safe Pack Pen Needles. To upgrade your pet's syringe or pen needle, simply print out this sheet and bring it to your vet.



Name	Needle Length	Gauge	Quantity	UPC	Built-in Sharps Container?
ULTICARE VETRX ULTIGUARD SAFE PACK U-40 INSULIN SYRINGES					
3/10 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-07263-4	Yes
3/10 mL/cc 12.7mm x 29G with Half Unit Markings	12.7mm (1/2")	29G (0.33mm)	100	3-57515-82003-7	Yes
1/2 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-07260-3	Yes
1/2 mL/cc 12.7mm x 29G with Half Unit Markings	12.7mm (1/2")	29G (0.33mm)	100	3-57515-82004-4	Yes
1 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-07261-0	Yes
ULTICARE VETRX ULTIGUARD SAFE PACK U-100 INSULIN SYRINGES					
3/10 mL/cc 8mm x 31G with Half Unit Markings	8mm (5/16")	31G (0.25mm)	100	3-57515-07436-2	Yes
3/10 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	3-57515-92001-0	Yes
3/10 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-07231-3	Yes
1/2 mL/cc 8mm x 31G	8mm (5/16")	31G (0.25mm)	100	3-57515-92002-7	Yes
1/2 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-07251-1	Yes
ULTICARE VETRX ULTIGUARD SAFE PACK PEN NEEDLES					
3/10 mL/cc 8mm x 31G with Half Unit Markings	8mm (5/16")	31G (0.25mm)	100	3-57515-07436-2	Yes

Insulin Syringe Ordering Information Continued:

UltiMed offers U-40 and U-100 insulin syringes in its UltiCare VetRx product lines.



Name	Needle Length	Gauge	Quantity	UPC	Built-in Sharps Container?
ULTICARE VETRX U-40 INSULIN SYRINGES					
3/10 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-09263-2	No
3/10 mL/cc 12.7mm x 29G with Half Unit Markings	12.7mm (1/2")	29G (0.33mm)	100	3-57515-82001-3	No
1/2 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-09260-1	No
1/2 mL/cc 12.7mm x 29G with Half Unit Markings	12.7mm (1/2")	29G (0.33mm)	100	3-57515-82002-0	No
1 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-09261-8	No
ULTICARE VETRX U-100 INSULIN SYRINGES					
3/10 mL/cc 8mm x 31G with Half Unit Markings	8mm (5/16")	31G (0.25mm)	100	3-57515-09436-0	No
3/10 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-09230-4	No
1/2 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-09250-2	No
1/2 mL/cc 12.7mm x 28G	12.7mm (1/2")	28G (0.36mm)	100	3-57515-08250-3	No
1 mL/cc 12.7mm x 29G	12.7mm (1/2")	29G (0.33mm)	100	3-57515-09210-6	No
1 mL/cc 12.7mm x 28G	12.7mm (1/2")	28G (0.36mm)	100	3-57515-08210-7	No