

Safe

Safety

- » Needle safely contained behind shield, preventing accidental exposure.

Comfort of use

- » Silicone lubrication and thin wall technology for effortless glide and medication flow.



Needle locking system

- » Needle is automatically locked after use.

Simple

- **Easy attachment** » Universal fit: compatible with most diabetes pen injector devices.
- **Simple use** » Just twist on, inject and dispose.
- **Needle lengths** » Two needle sizes to meet your needs (5mm x 30G, 8mm x 30G).

UltiCare®

Safety Pen Needles

Product intended use:

- » UltiCare safety pen needles are sterile, single-use safety pen needles intended for use with insulin pens for the subcutaneous injection of insulin.
- » Safety shield automatically locks in place and avoids cross infection caused by accidental needlestick injuries from the patient's end of the needle.
- » For use by clinical and non-clinical users.

Available in two sizes:



UltiCare 5mm

5mm (3/16") x **30G** (0.3mm)
NDC/HRI# 08222-6100-58
Product# 61005
100 pen needles/box
10 boxes/case
1,000 pen needles/case



UltiCare 8mm

8mm (5/16") x **30G** (0.3mm)
NDC/HRI# 08222-6100-89
Product# 61008
100 pen needles/box
10 boxes/case
1,000 pen needles/case

Instructions for use

