



ULTIMED

INCORPORATED

**Manufacturer of premium quality
injection devices and related
medical products with a focus on
serving people and pets
with diabetes.**



A leader in safe sharps disposal.

For over 30 years, UltiMed Incorporated has manufactured premium quality needles that provide superior comfort to millions of people living with diabetes.

UltiMed is committed to improving the lives of our customers and communities. As a world leader in safe sharps disposal, we focus on the entire product life cycle through innovation of easy to use, all-in-one solutions.



UltiGuard Safe Pack Sharps Container & Mail-Back Disposal Kit

4

PEN NEEDLES & INSULIN SYRINGES

UltiGuard Safe Pack

6

PEN NEEDLES & INSULIN SYRINGES

UltiCare

8

PEN NEEDLES

UltiCare

10

INSULIN SYRINGES

UltiCare

12

SAFETY PEN NEEDLES & ALCOHOL SWABS

UltiCare

14

SAFETY SYRINGES

UltiGuard[®] SafePack

Sharps Container & Mail-Back Disposal Kit

Bilingual
packaging

The UltiGuard Safe Pack Sharps Container & Mail-Back Disposal Kit is the only complete & compliant 3-in-1 solution that meets sharps disposal regulations in all 50 states. The kit includes 100 pen needles or insulin syringes, an FDA-cleared sharps container, and a prepaid return shipping box for shipment to a medical waste treatment facility for proper destruction. The UltiGuard Safe Pack Sharps Container & Mail-Back Disposal Kit requires no change in pharmacy workflow, and in most cases, patient copays will be the same as pen needles or syringes alone.








Pen Needles

100-COUNT

Micro: 4mm (5/32") × 32G	NDC# 08222-1254-39	
Mini: 5mm (3/16") × 31G	NDC# 08222-1255-38	
Mini: 6mm (1/4") × 32G	NDC# 08222-1256-20	
Short: 8mm (5/16") × 31G	NDC# 08222-1258-35	
Original: 12.7mm (1/2") × 29G	NDC# 08222-1251-25	

Insulin Syringes

100-COUNT

3/10 mL/cc 8mm (5/16") x 31G	NDC# 08222-1243-95	
1/2 mL/cc 8mm (5/16") x 31G	NDC# 08222-1245-93	
1 mL/cc 8mm (5/16") x 31G	NDC# 08222-1241-97	
3/10 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-1233-50	
1/2 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-1235-58	
1 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-1231-52	







UltiGuard® SafePack

Bilingual
packaging



The UltiGuard Safe Pack is the only all-in-one solution to provide premium, 100-count pen needles or insulin syringes and an FDA-cleared sharps container. In most cases, patient copays will be the same as pen needles or syringes alone.




Pen Needles

100-COUNT		
Micro: 4mm (5/32") × 32G	NDC# 08222-0954-35	
Mini: 5mm (3/16") × 31G	NDC# 08222-0955-34	
Mini: 6mm (1/4") × 31G	NDC# 08222-0956-33	
Mini: 6mm (1/4") × 32G	NDC# 08222-0956-26	
Short: 8mm (5/16") × 31G	NDC# 08222-0958-31	
Original: 12.7mm (1/2") × 29G	NDC# 08222-0951-21	







Insulin Syringes




100-COUNT		
3/10 mL/cc 6mm (1/4") x 31G	NDC# 08222-9101-03	
3/10 mL/cc 6mm (1/4") x 31G with half unit marking	NDC# 08222-9101-10	
1/2 mL/cc 6mm (1/4") x 31G	NDC# 08222-9101-27	
1 mL/cc 6mm (1/4") x 31G	NDC# 08222-9101-34	
3/10 mL/cc 8mm (5/16") x 31G	NDC# 08222-0743-93	
1/2 mL/cc 8mm (5/16") x 31G	NDC# 08222-0745-91	
1 mL/cc 8mm (5/16") x 31G	NDC# 08222-0741-95	
3/10 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-0733-58	
1/2 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-0735-56	
1 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-0731-50	

UltiCare pen needles are precision engineered to provide superior comfort. Every needle is individually inspected using laser technology, ensuring a more comfortable injection.



Pen Needles

100-COUNT		
Micro: 4mm (5/32") × 32G	NDC# 08222-7100-31	
Mini: 5mm (3/16") × 31G	NDC# 08222-7101-30	
Mini: 6mm (1/4") × 31G	NDC# 08222-7100-48	
Mini: 6mm (1/4") × 32G	NDC# 08222-7101-47	
Short: 8mm (5/16") × 31G	NDC# 08222-7100-55	
Original: 12.7mm (1/2") × 29G	NDC# 08222-7100-62	



















	30-COUNT	50-COUNT	
Micro: 4mm (5/32") × 32G	NDC# 08222-7100-79	NDC# 08222-0954-59	
Mini: 6mm (1/4") × 31G	NDC# 08222-7100-86	NDC# 08222-0956-57	
Short: 8mm (5/16") × 31G	NDC# 08222-7100-93	NDC# 08222-0958-55	

UltiCare insulin syringes combine proven comfort and accuracy with features like triple bevel silicone lubricated needle for effortless glide, bold barrel markings for accurate dosing, and large flanges for better grip and control.



Insulin Syringes

100-COUNT

3/10 mL/cc 6mm (1/4") x 31G	NDC# 08222-9100-66	
3/10 mL/cc 6mm (1/4") x 31G with half unit marking	NDC# 08222-9100-73	
1/2 mL/cc 6mm (1/4") x 31G	NDC# 08222-9100-80	
1 mL/cc 6mm (1/4") x 31G	NDC# 08222-9100-97	
3/10 mL/cc 8mm (5/16") x 31G	NDC# 08222-0943-91	
1/2 mL/cc 8mm (5/16") x 31G	NDC# 08222-0945-99	
1 mL/cc 8mm (5/16") x 31G	NDC# 08222-0941-93	
3/10 mL/cc 8mm (5/16") x 30G	NDC# 08222-0933-94	
1/2 mL/cc 8mm (5/16") x 30G	NDC# 08222-0935-92	
1 mL/cc 8mm (5/16") x 30G	NDC# 08222-0931-96	
3/10 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-0933-56	
1/2 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-0935-54	
1 mL/cc 12.7mm (1/2") x 30G	NDC# 08222-0931-58	
3/10 mL/cc 12.7mm (1/2") x 29G	NDC# 08222-0923-97	
1/2 mL/cc 12.7mm (1/2") x 29G	NDC# 08222-0925-95	
1 mL/cc 12.7mm (1/2") x 29G	NDC# 08222-0921-99	
1/2 mL/cc 12.7mm (1/2") x 28G	NDC# 08222-0825-89	
1mL/cc 12.7mm (1/2") x 28G	NDC# 08222-0821-83	





UltiCare safety pen needles are precision engineered with advanced technology to help keep clinicians safe while providing patients with a more comfortable injection. Compatible with most pen injector devices in the United States and tested in accordance with ISO standards for compatibility.


UltiCare alcohol swabs are individual pads saturated in 70% isopropyl alcohol for optimum anti-bacterial action.



Safety Pen Needles

100-COUNT		
Mini: 5mm (3/16") × 30G	NDC# 08222-6100-58	
Short: 8mm (5/16") × 30G	NDC# 08222-6100-89	

Alcohol Swabs



100-COUNT	200-COUNT
NDC# 08222-4400-06	NDC# 08222-4400-13 

UltiCare safety syringes are designed to protect patients and healthcare professionals from accidental needle stick injuries. The insulin and tuberculin safety syringes with attached needle provide low dead space design to minimize medication waste.




Insulin Safety Syringes (WITH ATTACHED NEEDLE)

100-COUNT

1/2 mL/cc 12.7mm (1/2") x 29G	NDC# 08222-0325-91	
1 mL/cc 12.7mm (1/2") x 29G	NDC# 08222-0321-95	








Tuberculin Safety Syringes (WITH ATTACHED NEEDLE)

100-COUNT

1 mL/cc 12.7mm (1/2") x 28G	NDC# 08222-6300-25	
1 mL/cc 12.7mm (1/2") x 27G	NDC# 08222-6300-18	
1 mL/cc 15.9mm (5/8") x 27G	NDC# 08222-6300-01	
1 mL/cc 15.9mm (5/8") x 25G	NDC# 08222-2515-89	
1 mL/cc 25.4mm (1") x 25G	NDC# 08222-2511-07	

Luer Lock Safety Syringes (WITH DETACHABLE NEEDLE)

100-COUNT

3 mL/cc 25.4mm (1") x 25G	NDC# 08222-6300-87	
3 mL/cc 15.9mm (5/8") x 25G	NDC# 08222-6300-70	
3 mL/cc 25.4mm (1") x 23G	NDC# 08222-6300-63	
3 mL/cc 38.1mm (1-1/2") x 22G	NDC# 08222-6300-56	
3 mL/cc 25.4mm (1") x 22G	NDC# 08222-6300-49	
3 mL/cc 38.1mm (1-1/2") x 21G	NDC# 08222-6300-32	
3 mL/cc Syringe Only	NDC# 08222-6301-48	



ULTIMED
INCORPORATED

ultimedinc.com · 1-877-854-3434